


DogLeggs Carpal Support is easiest to apply and adjust when dog is lying down and leg and paw are in a neutral position. The DogLeggs Carpal Support can be used on either leg. We suggest applying Carpal Support with the Tabs wrapping around the leg and finishing toward the outside of the leg (see example).

*\*Dogs with auto-immune issues or long term prednisone therapy, must wear a breathable, lined Carpal Support .*

Care: Machine wash gentle, air dry

### Application of Carpal Support



2104 Thomas View Road, Reston VA 20191

800-313-1218 [www.dogleggs.com](http://www.dogleggs.com) fax: 703-391-9333