

Surgi-Sox™ Animal Bandaging System



The DogLeggs Surgi-Sox™ Animal Bandaging System provides coverage, protection and compression for wounds, abrasions, skin lesions and incisions located on the torso, forelimbs, head and ears without restriction of movement. Quick and easy to slip on, these products can often replace traditional bandaging methods and the Elizabethan collar.

All components of the Surgi-Sox™ product line are made from a proprietary material. This 4-way stretch textile is highly breathable, wicks away moisture, is plated to prevent the fabric from sticking to wounds and has anti-bacterial and anti-microbial properties woven within. These properties withstand hundreds of washings. VELCRO® Brand fasteners secure the bandage. The product is machine washable.

Surgi-Sox™ Aural Compression Bandage

This revolutionary alternative for head and ear bandaging uses a compressive quality that prevents damage to the ear during head shaking. The Surgi-Sox™ Aural Compression Bandage slips on and VELCRO® Brand fasteners provide secure adjustments just in front and behind the ears.

Indications for Use:

- 🐾 Ear Hematomas
- 🐾 Wounds, abrasions, incisions on or around the ear
- 🐾 Total Ear Canal Ablation
- 🐾 Shaking Splatter following ear medication application

**Surgi-Sox™
Torso Bandage**
(See reverse side for
product details)



**Surgi-Sox™
Aural Compression
Bandage**
(See reverse side for product
measurements)

Sizing for Aural Compression Bandage:



1 Measure around head just in front of ears.

Aural Compression Bandage		
Size	Inches	SKU#
S	13"-16"	6500.03
M	15"-18"	6500.05
L	17"-21"	6500.07
XL	20"-25"	6500.08



Surgi-Sox™ Torso Compression Bandage

The Surgi-Sox™ Torso provides coverage, protection and compression of the torso. Compression can help prevent seroma formation following surgery and can also benefit animals with certain behavioral issues. Use of the Surgi-Sox™ Torso can often replace traditional bandaging methods and the Elizabethan Collar. Surgi-Sox™ products are easy to apply. Improved patient comfort reduces licking, chewing and scratching at wounds. The simplicity of the product results in greater compliance. The level of compression can be adjusted based on the size selected. Our sizing chart is based on a mild to moderate compression and can be used as a guideline if more or less compression is desired.

Sizing for Torso Compression Bandage:



1 Measure the circumference of the chest immediately behind the front legs (at its deepest point).

		Surgi-Sox™ Torso	Surgi-Sox™ Torso for the Forelimb Amputee	Surgi-Sox™ Torso with VELCRO® (down the back)	Surgi-Sox™ Torso with Leggings (front)	Surgi-Sox™ Torso with Leggings (front) and VELCRO® (down the back)
Size	Inches	SKU#	SKU#	SKU#	SKU#	SKU#
XS	8.5"-10.5"	6000.02	6100.02	6200.02	6300.02	6400.02
S	10"-12.5"	6000.03	6100.03	6200.03	6300.03	6400.03
S/M	11.5"-14"	6000.04	6100.04	6200.04	6300.04	6400.04
M	14"-17.5"	6000.05	6100.05	6200.05	6300.05	6400.05
M/L	18"-22.5"	6000.06	6100.06	6200.06	6300.06	6400.06
L	24"-30"	6000.07	6100.07	6200.07	6300.07	6400.07
XL	26"-32.5"	6000.08	6100.08	6200.08	6300.08	6400.08
XXL	28"-35"	6000.09	6100.09	6200.09	6300.09	6400.09
3XL	31"-38.5"	6000.10	6100.10	6200.10	6300.10	6400.10

Indications for Use:

- Wounds, abrasions, incisions
- Skin lesions associated with dermatological disease
- Chronic obsessive licking of the torso
- Surgically placed drain tubes/ Jackson Pratt drains
- Certain behavioral issues (anxiety/stress)



Surgi-Sox™ Torso Compression Bandage



Surgi-Sox™ Torso for the Forelimb Amputee



Surgi-Sox™ Torso with VELCRO® (down the back)



Surgi-Sox™ Torso with Leggings (front)



Surgi-Sox™ Torso with Leggings (front) and VELCRO® (down the back)