

VEST System






The DogLeggs VEST is the platform for the DogLeggs Ehmer Sling and Cardiac Monitoring Solution.

Ehmer Sling

The VEST with Ehmer Sling is a superior alternative to bandaging. The products are machine washable and reusable for multiple patients. The VEST with Ehmer Sling eliminates the trauma to the skin caused by traditional bandaging, improving patient comfort and recovery. The system is easy to apply and can be used for either hind limb. Rehabilitation therapy can be performed with ease throughout recovery.


Indications for Use:

-  Post Closed/Open Reduction of Craniodorsal Coxofemoral Luxations
-  Post Operative Management of Acetabular, Femoral Head and Neck and certain Femoral Fractures
-  To Encourage Weight-bearing in the Contralateral Limb for Rehabilitative Therapy

Cardiac Monitoring

The VEST for Cardiac Monitoring consists of the DogLeggs VEST, Device Case and Additional Stabilization Strap. The VEST has a small opening along the dorsal aspect on the right side to connect the leads from the animal to the device. The Device Case secures to the VEST using VELCRO® Brand fasteners and the Additional Stabilization Strap provides further securement of the Device Case to the VEST. This system protects the device and the VEST protects the leads. The Device Case is made to fit the specific device model. One Device Case may be utilized for multiple VESTs.

Indications for Use:

-  Need for consistent reliable monitoring of animals in the clinic or at home



VEST for Ehmer Sling
(See reverse side for product measurements)




VEST for Cardiac Monitoring
(See reverse side for product measurements)

Specifications:

The DogLeggs VEST and supplemental components are constructed from a lightweight material that is VELCRO® Brand fastener compatible. This material is highly breathable and has wicking properties. The material does not pull the fur and is impregnated with a natural fiber that reduces odor and inhibits bacterial and fungal growth. VELCRO® Brand fasteners are used for closure and attachments. The product and supplemental components are machine washable.

VEST Sizing for Ehmer Sling and Cardiac Monitoring:

- 1 Measure the circumference of the chest immediately behind the front legs (at its deepest point).



Size	Inches	VEST for Cardiac Monitoring	VEST with Ehmer Sling
		SKU#	SKU#
XXXS	10"-13"	3200.00	3300.00
XXS	14"-16"	3200.01	3300.01
XS	17"-19"	3200.02	3300.02
S	20"-22"	3200.03	3300.03
S/M	23"-25"	3200.04	3300.04
M	25"-27"	3200.05	3300.05
M/L	28"-34"	3200.06	3300.06
L	35"-39"	3200.07	3300.07
XL	40"-43"	3200.08	3300.08
XXL	44"-58"	3200.09	3300.09

Cardiac Monitoring Device Case
SKU#
5100.99

